



DNA LYSIS BUFFER FROM WHOLE BLOOD

User Guide

Catalog Number	Pack Size
ZT-BLB-100	100 ML
ZT-BLB-500	500 ML

Content and Storage Condition

Content	Shipping Condition	Storage Condition
DNA LYSIS BUFFER FROM WHOLE BLOOD	RT	RT

Product Description

The DNA Lysis Buffer From Whole Blood provides an important protection against the high nuclease activity found in blood samples, while ensuring the rapid breakdown of hemoglobin and other protein components.

Concentration

1X

Appearance

Liquid

Warning

Do not preserve the product when the package is damaged.

